

# Disposable Blood Filter

## Indications

suitable for filtering whole blood, red blood cells, platelets, concentrated granulocytes, plasma, condensed proteins, Viaspan refrigeration fluid

## Specifications

Code	Description	Unit
BFA40-S	Filter 40µm without infusion line	10/pack
BFA20-S	Filter 20µm without infusion line	10/pack
BFA40-L	Filter 40µm with infusion line	10/pack
BFA20-L	Filter 20µm with infusion line	10/pack

## Features

- ※ Fast filtration, can filter red blood cells up to 2000ml
- ※ Easy to defoaming
- ※ Easy to observe
- ※ High filtration efficiency
  - Can filter out more than 95% of impurities above 20 microns (BFA20 type)
  - Can filter out more than 95% of impurities above 40 microns (BFA40 type)
- ※ Maximum working pressure: 400mmHg
- ※ Polyester monofilament filter element

